


Cardiac MRI




O. Wieben
Depts. of Radiology and Medical Physics
University of Wisconsin-Madison

Cardiac MRI



Cardiac Function Velocity Mapping Tagging Valve Imaging
 Coronary Artery Imaging Perfusion Viability

What is special about cardiac MR?



Dealing with motion

- Cardiac motion
- Breathing

Adjust sequences accordingly

- ECG Gating
- Breathhold – Navigator
- Rapid Imaging

X-ray DSA Matrix: 512x512, 30 frames/s

Outline

Physiological Monitoring
 ECG Gating
 Breathholds and Navigators


Sequences
 Bright Blood: Gradient echo sequences
 Phase Contrast Imaging
 Dark Blood: Double IR
 Perfusion
 Viability
 Tagging
 Coronary Imaging

1.5T vs. 3T

Not covered: CE-MRA and NCE-MRA

Cardiac MRI

Cardiac Function, Anatomy	Clinically applied
Myocardial Perfusion	
Myocardial Viability, Infarct Size	Getting there ?? 64 Slice CT competition ?
Valvular Heart Disease	
Tagging	Predominantly Research
Coronary Angiography	
Plaque Imaging and Composition	
T1 mapping	
Interventional Imaging	
Stem Cell Tracking	



ELSEVIER European Journal of Radiology 65 (2006) 15–28 **EJR**
 www.elsevier.com/locate/ejrad

Review
Cardiac MRI of ischemic heart disease at 3 T: Potential and challenges
 Oliver Wieben^{a,b,*}, Christopher Francois^{a,1}, Scott B. Reeder^{a,b,c,d,2}


^a Department of Radiology, University of Wisconsin, Madison, WI 53702-3252, United States
^b Department of Medical Physics, University of Wisconsin, Madison, WI 53702-3252, United States
^c Department of Biomedical Engineering, University of Wisconsin, Madison, WI 53702-3252, United States
^d Department of Medicine, University of Wisconsin, Madison, WI 53702-3252, United States

Received 30 October 2007; accepted 30 October 2007

Getting ready



Position coil correctly
Get a good ECG signal

- Vector ECG
- Wireless ECG



Physiological Monitoring

Breathing Motion



Temporal resolution: 0.3s
Playback: 3x accelerated


Deal with Breathing Motion

Breathhold



- End-inspiration
- End-expiration
- Limits scan time to 10-20s
- Not feasible for some patients

Free breathing

- Bellows
- Navigators
- Beat the motion – realtime MRI




Navigator Position



Stuber, Manning et al. JMRI 10:790-799 (1999)

Navigator (Bellow) Analysis



Gating window: 4-6 mm
Widen window with slice tracking
Efficiency for coronaries: 40%

Navigator – Potential Issues


Low temporal resolution

- Once per RR interval

Interrupt steady state sequences

- Do not directly measure motion of the heart
 - Measure diaphragm position instead

Limited precision



Cine vs. Non-cine


Cine

Cardiac function
Phase contrast
Tagging



Non-cine

Coronaries
Viability
Double & Triple IR




Shortest Scan duration



Example

Heartrate = 60bpm
RR cycle = 1s
TR = 3.3ms
Acquisition matrix: 256x192
Parallel Imaging: acceleration x 2

Shortest Scan duration




Example

Heartrate = 60bpm
RR cycle = 1s
TR = 3.3ms
Acquisition matrix: 256x192
Parallel Imaging: acceleration x 2

Total scan duration:
96*TR = 310 ms
Temporal resolution: 310 ms


Shortest Scan duration



Example

Heartrate = 60bpm
RR cycle = 1s
TR = 3.3ms
Acquisition matrix: 256x192
Parallel Imaging: acceleration x 2

Total scan duration:
96*TR = 310 ms
Temporal resolution: 310 ms




Myxoma: Fast localization and Acquisition

Real-time Parameters Scan and Save Parameters















Real Time series



"Scan & Save" series

Image courtesy Dr Sablayrolles, Centre Cardiologique du Nord





Off-resonance




Off-Resonance Artifact




Corrected by re-shimming

Banding Artifacts with bSSFP




avoid bandartifacts!
TR < 5 ms (1.5 T)

Effect of TR




Adjusting Center Frequency



CF = -50




Center frequency = 0










CF = +70

Standard CMR imaging planes




LV function







Dark Blood Imaging



Double IR (Blood Suppression)




Double IR: Short and Long Axis Views



Double IR: Vessel Wall Imaging



Double IR: Valve Imaging




Summary: Dark Blood Cardiac Sequences

- Employ ECG-gated, segmented, fast spin echo sequences
- Allow acquisition of multiple slices and a single cardiac phase for each slice
- May employ double inversion for blood suppression or triple inversion for blood and fat suppression
- Provide better soft tissue contrast than GRE sequences
- Used for assessing anatomy/cardiac structure – masses, thrombi, acute or remote myocardial infarctions, valves, vessel walls, etc.


Myocardial Perfusion

Myocardial Perfusion




- First pass imaging during injection of contrast agent
- Sequences:
 - multiecho GRE,
 - single echo GRE with parallel imaging
- Inject 0.05 mmol/kg at 4-5 ml/second
- 5-6 Slices in SA
- Stress/Rest using adenosine (dipyrimadole) and reversing with aminophylline
- T_2^* artifact especially in septum

Stress Myocardial Perfusion




Nuclear Medicine




Rest MRI

Stress MRI

Saturation Recovery Sequence




1 vs 4 echo train length (ETL)




Perfusion - Artifacts


- Clinical routine at 1.5 T:
 - Patient is imaged under stress, then reverted and imaged again
 - FGRET sequence



1st Scan: stress FGRET



2nd Scan: rest FGRET



Viability Imaging
Late Gd Enhancement (LGE)
Myocardial Delayed Enhancement (MDE)

