

Non-Cartesian MRI

Oliver Wieben, Ph.D.
 Depts. of Medical Physics & Radiology
 University of Wisconsin - Madison

Image Encoding – The early days

First 'image': Lauterbur, 16 Mar 1973, Nature
 2 water-filled capillary tubes

First human body MRI Chest:
 Damadian et al. 1977

First human body MRI
 with line-scan:
 Mansfield et al, 1978

First human MRI: finger:
 Mansfield et al, 1976

Spinwarp Imaging

Figure 1. Spin warp imaging pulse sequence.

W. A. Edelstein et al.
 Spin warp NMR imaging and applications to human whole-body imaging
 Physics in Medicine and Biology 25: 751-6, 1980

Spinwarp Imaging







W. A. Edelstein et al., 1980
 0.04 T; Image matrix: 64x64
 Scan time: 128 s





Data Acquisition – k-space

2D Gradient Echo

Data Acquisition – k-space

2D Gradient Echo









Alternative Trajectories

Radial Sampling
Every k-space line traverses through center of k-space

Propeller
Thin k-space slabs rotating around origin

Spiral imaging
Long readouts in spiral trajectory

More exotic trajectories
Cones, Rosettes, Shells, Hybrids, ...

Design Goals

Reduced total scan times / larger coverage
Improved efficiency (higher readout/TR ratio)
Undersampling with tolerable artefacts

Short TE (ZTE, UTE)
Reduced flow artefacts (small 1st moment around TE)
Reduced intravoxel dephasing and displacement

Dealing with Motion
Reduced motion artefacts, Self gating, Double Gating

'Random' sampling
Incoherent sampling artefacts for compressed sensing





Alternative Trajectories

Radial Sampling
Every k-space line traverses through center of k-space

Propeller
Thin k-space slabs rotating around origin


Spiral imaging
Long readouts in spiral trajectory

More exotic trajectories

Fast Three-Dimensional k-Space Trajectory Design Using Missile Guidance Ideas

Roberto Mir, Andres Gussalaga, Juan Spiniak, Marcelo Guarini, and Pablo Irarrazaval*




Alternative Trajectories


Radial Sampling
Every k-space line traverses through center of k-space

Propeller
Thin k-space slabs rotating around origin

Spiral imaging
Long readouts in spiral trajectory

More exotic trajectories



Design goals

Design Limitations





Alternative Trajectories

Radial Sampling
Every k-space line traverses through center of k-space

Propeller
Thin k-space slabs rotating around origin

Spiral imaging
Long readouts in spiral trajectory


More exotic trajectories

Design goals

Design Limitations

Spiral Imaging






FIG. 2. Fast interleaved spiral pulse sequence. There is a spectral-spatial pulse for fat suppression followed by spiral readout gradients. The sequence is designed to image a 20 cm FOV in 20 excitations, with a resolution of 1.08 mm by 1.08 mm (156 × 156 pixels). The readout gradients corresponding to one interleaf are shown; the other interleaves are generated by rotating the gradients with a hardware rotation board.


Prior GE System MR 4, 1.5T: whole body (48 cm) zoom (35 cm)

max gradient strength = 23 mT/m 40 mT/m


max gradient slew rate = 77 mT/m/s 150 mT/m/s







C. Meyer et al., MRM 1992







2D Radial Trajectory



$P(k_r, \theta)$







Radial Acquisitions

Rapid imaging with undersampling


- Interventional MR *K. Scheffler, MRM 1998*
- Contrast-enhanced MRA *S. Weiss, MRI 1999*
- bSSFP cardiac cine *D. Peters, MRM 2000*
- A. Larson, MRM 2001*
- T. Schaeffter, MRM 2001*
- Shankaranarayanan, Radiology 2001*

Robust to motion artifacts

- Motion is 'averaged' *G. Glover, MRM 1992*
- No ghosting artifacts
- Navigational information *N. Gai, Med Phys 1996*

UTE (Ultrashort TE imaging)

- Minimal TE times for short T2 *G. Glover, JMRI 1992*
- CSI imaging *J. Du, MRM 2007*



Non-Cartesian – 3D Trajectories




FIG. 1. 3DFT trajectory.




FIG. 2. Cylindrical 3DFT trajectory.

<p>10,000 excitations (8 min 20 s)</p> <p>Prep time = 1.43 ms</p> <p>TE = 2.3 ms</p> <p>Significant first moment at TE</p> <p>Flow Comp would be best</p>	<p>7,854 excitations (6 min 33 s)</p> <p>Prep time = 1.43 ms</p> <p>TE = 2.3 ms</p> <p>Significant first moment at TE</p> <p>Flow Comp would be best</p>
---	--

100x100x100 pixels
TR = 50 ms
Readout time = 14.3 ms

Max $M_z / 2s \leq 1$ G/cm
Max $G/cm/ms$

From P. Irarrazabal and D. Nishimura, MRM 1995

Non-Cartesian – 3D Trajectories




FIG. 5. Stack of spiral trajectory.

900 excitations (45 s)
Prep time = 1.43 ms
TE = 0 ms
Low first moment at TE
but slice encoding




FIG. 6. Spherical stack of spiral trajectory.

628 excitations (31 s)
Prep time = 1.43 ms
TE = 0 ms
Low first moment at TE
but slice encoding

From P. Irarrazabal and D. Nishimura, MRM 1995

Non-Cartesian – 3D Trajectories




FIG. 7. Shells trajectory.

629 excitations (31 s)
Prep time = 1.43 ms
TE = 0 ms
Low first moment at TE
but 'slice' encoding (pole)




FIG. 8. 3D projection-reconstruction trajectory.

31,416 excitations (26 min 10 s)
Prep time = 0 ms
TE = 0 ms
zero first moment at TE
excellent flow properties

From P. Irarrazabal and D. Nishimura, MRM 1995

Non-Cartesian – 3D Trajectories




FIG. 10. Yam trajectory.




FIG. 9. Spiral-PR hybrid trajectory.






FIG. 11. Cones trajectory.

From P. Irarrazabal and D. Nishimura, MRM 1995


Radial Sampling



2D




3D hybrid – Stack of stars



3D radial - VIPR

VIPR

Vastly Undersampled Isotropic PROjection Imaging



3D Radial acquisition

Nyquist: $N_p = \pi/2 N_{\text{readout}}^2$


Large FOV

Isotropic spatial resolution


Undersampling

A. Barger et. al, MRM 2002


VIPR



Vastly undersampled
Isotropic Projection
Reconstruction



A. Barger et. al, MRM 2002



PC VIPR – Cranial

Normal Volunteer

PC VIPR Parameters

- Dual Echo
- FOV: 20 x 20 x 20 cm
- Res: 0.6 x 0.6 x 0.6 mm
- 9000 Projections (36x)
- TR=15.9
- Bandwidth = 31.25
- VENC = 50 cm/s
- 5:07 min Scan Time

Same Cartesian PC

- 48+ min Exam (Partial)

Same TOF

- 24+ min Exam (Partial)




Image Reconstruction





Image Reconstruction

Projection Reconstruction

- Matrix inversion
- Iterative Reconstruction
- Filtered back projection (Radon)
 - (magnitude or complex)









Godfrey Hounsfield

Gridding (Regridding)

- Weight data according to sampling density
- grid polar (spherical) samples onto Cartesian grid
- 2D or 3D inverse Fourier Transform
- Divide by FT of the interpolation kernel to remove image apodization from convolution


- Non-uniform FFT (NUFFT)



Gridding

1. Weight data according to sampling density
As simple as p for 2D or p^2 for 3D
2. Grid onto Cartesian Coordinates
Kaiser-Bessel or 2D/3D separable triangular convolution kernel
3. 2D or 3D inverse Fourier Transform
4. Divide by FT of the interpolation kernel to remove image apodization from convolution

-> much faster than Radon in 3D



J. Jackson et al., IEEE TMI 1991
M. Lauzon and B. Rutt, MRM 1998
J. Pipe, MRM 2000

Ultra-Short TE (UTE) Imaging

Tissues with short T2

T2 of protons in water
Usually clinical range


T2 of protons in proteins and in water tightly bound to protein
10 μ s range

Standard MR sequences
Detectable T2s: 10 ms, modern systems: 2 ms


Table 2
Tissues With a Majority of Short T₂ Components

Meninges (dura)	Falx	Tentorium
Membranes	Capsules	Bands
Retinaculi	Septae	Fasciae
Sheaths	Nails	Hair
Aponeuroses	Tendons	Ligaments
Menisci	Labrii	Periosteum
Bone	Dentine	Enamel


2D Radial Trajectory




$P(k_r, \theta)$



2D Radial Trajectory




Full echo radial




Fractional echo radial







2D Radial Trajectory









Full echo radial



UTE radial





Trajectory Errors

The desired and the actual trajectory differ

Gradient anisotropy

- different system delays for x, y, and z gradients
- > echo shifts: first order trajectory errors


B. Aldefeld, MRM 1998
S. Reeder, MRM 1999

Eddy currents

Short TRs and steep gradients


- > pre-emphasis corrections are for longer time constants
- > first and higher order trajectory errors

Data Acquisition - 2D GRE



Off-center echoes

- Leads to phase twist in image space
- FT imaging: magnitude is not affected
- Problems in oblique EPI and non-cartesian sampling



Trajectory Errors

The desired and the actual trajectory differ

Gradient anisotropy

- different system delays for x, y, and z gradients
- > echo shifts: first order trajectory errors

B. Aldefeld, MRM 1998
S. Reeder, MRM 1999

Eddy currents

Short TRs and steep gradients

- > pre-emphasis corrections are for longer time constants
- > first and higher order trajectory errors

Solutions

Correct first order trajectory errors

- Measure and analyze antiparallel echoes
- Correct

D. Peters, ISMRM 2003
O. Wieben, ISMRM 2003


During gridding

- During acquisition: modify gradient area


Measure actual trajectory

Duyn, JMR 1998
B. Dale, ISMRM 2003


Results: Ball Phantom



Results: High Resolution Phantom





K-space Trajectory Measurement



- Slice thickness: 1mm~1.5 mm
- Off-center distance (D_r): 3 cm

- Can be performed during the imaging scan

Duyn, et al. JMR 132 150, 1998 B. Dale, 10th ISMRM 2003, 2334



Non-Cartesian PC MRI

Benefits

- + Accelerated Imaging
- + High spatial resolution (isotropic)
- + Favorable Motion artefacts
- + UTE Flow Imaging
- + (Large Coverage)

Non-standard Acquisition

Reconstruction is more complex

- Gridding
- → Longer recon times
- Parallel Imaging Recon is challenging

Sensitivity to Artefacts

- Trajectory errors (calibration required)
- Off-resonance

Artefacts

- Streaks/Blurring unfamiliar to clinicians
- Displacement in readout (changing dir)

Summary

Cartesian

- Robust image quality
- Proven
- Fast reconstruction
 - works well with parallel imaging
- Tailor Image Volume

Non-Cartesian

- Can overcome some of current limitations
 - Scan time / spatial resolution / coverage
 - Short TE
 - 'Random Sampling'
- Acquisition more complex
 - calibrations
- Recon more complex
 - Feasible (see Propeller)