






MP 710 – Lecture 3 Sept. 2017



O. Wieben, Ph.D.
Depts. of Medical Physics & Radiology
University of Wisconsin - Madison





MRI – Selected Online Resources

<http://www.cis.rit.edu/htbooks/mri/> <https://www.mriquestions.com/index.html>

MR Signal – Quantum Mechanics



Is Quantum Mechanics Necessary for Understanding Magnetic Resonance?

LARS G. MANDON
David Rossier Center for Magnetic Resonance, Cambridge University Hospital, Madbury, Chesham

ABSTRACT: Educational material introducing magnetic resonance (MR) typically contains sections on the underlying principles. Unfortunately, the explanations given are often unnecessarily complicated or even wrong. MR is often presented as a phenomenon that necessitates a quantum mechanical explanation whereas in reality it is a classical effect. In a commentary of the common sense approach to classical mechanics, this might be not true, but these books have an attempt to challenge common misperceptions. There are sections and chapters on fundamental concepts like spin, as a result, some students' first encounters with MR are often obscured by explanations that make the subject difficult to understand. Typical students are addressed and alternative theories explanations are provided. © 2006 Wiley Periodicals, Inc. Concepts Magn Reson Part 1: Clin 10: 405-406, 2006

Nuclei with spin – relevant to MR

The shell model for the nucleus tells us that nucleons, just like electrons, fill orbitals. When the number of protons or neutrons equals 2, 8, 20, 28, 50, 82, and 126, orbitals are filled. Because nucleons have spin, just like electrons do, their spin can pair up when the orbitals are being filled and cancel out. Almost every element in the periodic table has an isotope with a non zero nuclear spin. NMR can only be performed on isotopes whose natural abundance is high enough to be detected, however some of the nuclei which are of interest in MRI are listed below.

Nuclei	Unpaired Protons	Unpaired Neutrons	Net Spin	γ (MHz/T)
^1H	1	0	1/2	42.58
^2H	1	1	1	6.54
^{13}C	1	0	1/2	17.25
^{23}Na	1	2	3/2	11.27
^{14}N	1	1	1	3.08
^{31}P	0	1	1/2	10.71
^{19}F	1	0	1/2	40.08

Source: <http://www.cis.rit.edu/htbooks/mri/>

31P, 13C, 14N, spectroscopy
13C hyperpolarized gas
23 Na – sodium imaging – 90% extracellular

Abundance

Natural abundance of an isotope is the fraction of nuclei having a given number of protons and neutrons, or atomic weight

Element	Symbol	Natural Abundance
Hydrogen	^1H	99.985
	^2H	0.015
Carbon	^{13}C	1.11
	^{14}N	99.63
Nitrogen	^{15}N	0.37
Sodium	^{23}Na	100
Phosphorus	^{31}P	100
Potassium	^{39}K	93.1
Calcium	^{40}Ca	0.145

Source: <http://www.cis.rit.edu/htbooks/mri/>


Biological abundance - fraction of one type of atom in the human body

Element	Biological Abundance*
Hydrogen (H)	0.63
Sodium (Na)	0.00041
Phosphorus (P)	0.0024
Carbon (C)	0.094
Oxygen (O)	0.26
Calcium (Ca)	0.0022
Nitrogen (N)	0.015

* Calculated from


Source: <http://www.cis.rit.edu/htbooks/mri/>

Magnetic moment: w and w/u external field




Source: <http://www.mhhe.com/physsci/chemistry/carey/student/olc/ch13nmr.html>

Magnetic moment: w and w/o external field




Source: <http://www.mhhe.com/physsci/chemistry/carey/student/olc/ch13nmr.html>

Magnetic moment: w and w/o external field



Source: <http://www.mhhe.com/physsci/chemistry/carey/student/olc/ch13nmr.html>

Nuclear Spin



Source: <https://www.quora.com/How-are-atomic-nuclei-and-magnets-similar>

Magnetic Moment in external field





Figure 3.2 Nuclear magnetic moment vectors (a) pointing in random directions and (b) aligned in the direction of an external magnetic field.

Right Hand Rule – Cross Product



R. Hewitt - Original publication: ManufacturingET.org Immediate source: <http://www.manufacturinget.org/wp-content/uploads/2011/12/right-hand-rule.png>

Magnetic force on a moving charged particle
The direction of the cross product may be found by application of the right hand rule as follows: Using your right hand, Point your index finger in the direction of the first vector A. Point your middle finger in the direction of the second vector B. Your thumb will point in the direction of the cross product C.

From http://en.wikipedia.org/wiki/Right-hand_rule

Precession





Figure 3.3 Precession of a nuclear spin about an external magnetic field is similar to the wobbling of a spinning top in a gravitational field.

Source: https://en.wikibooks.org/wiki/Basic_Physics_of_Nuclear_Medicine/MRI_%20Nuclear_Medicine

'Left Hand Rule'




Right Hand: Prediction of direction of field (B), given that the current I flows in the direction of the thumb
From http://en.wikipedia.org/wiki/Right-hand_rule

which is known as the Larmor frequency. Second, precession of \vec{J} about \vec{B}_0 is clockwise if observed against the direction of the magnetic field. In practice, it is easy to determine the precession direction using the **left-hand rule**. That is, if the left thumb points in the direction of \vec{B}_0 , nuclear precession follows the direction of other fingers.
In addition to Eq. (3.16), nuclear precession can be described by an angular velocity vector defined as


$$\vec{\omega} = -\gamma \vec{B}_0 = -\omega_0 \hat{k}$$

Diagram illustrating a quantum circuit for a controlled NOT gate. The circuit starts with a qubit in the B_0 off state. It then splits into two paths, each with a different energy level: $E = \frac{1}{2} \hbar B_0$ (upper path) and $E = -\frac{1}{2} \hbar B_0$ (lower path). The paths recombine at a point labeled B_0 on.




3

RF system – single TR coil



from McFall, 1997

RF system – separate TR coils



from McFall, 1997

Birdcage Coil



MR Signal Detection


M: Magnetization

$$\vec{M} = M_0 \vec{e}_z$$

f_0 : Larmor Frequency

$$f_0 = \frac{\gamma}{2\pi} B_0$$

$$= 63.86 \text{ MHz}$$




B_1 = Transverse Field


$$B_1(t) = \hat{B}_1 \sin(2\pi f_0 t)$$







$$\hat{B}_1 = 20 \text{ } \mu\text{T} = 0.2 \text{ G}$$


Birdcage Coil – RF Excitation



Linear birdcage volume coil







Total Digital Imaging (TDI)

The SIGNA™ Architect offers startling advances in imaging and a total imaging war with TDI.

GE's Direct Digital Interface (DDI) employs an independent analog-to-digital converter to digitize signals from each of 128 RF channels, eliminating unnecessary noise enhancement.

Digital Micro Switching (DMS) technology represents a revolutionary advance in RF coil design by replacing analog blocking circuits with intelligent Micro Electro-Mechanical Switches (MEMS).

The SIGNA™ Architect comes prepared for **Digital Surround Technology (DST)**. DST combines signals from every coil element. The exceptional SNR and sensitivity of the high-density surface coils are combined with the superior homogeneity and deeper signal penetration of the integrated RF body coil.

http://www3.gehealthcare.com/en/products/categories/magnetic_resonance_imaging/3-0/signa_architect


GE - Air Coils

AIR Technology


- 80% more flexible design
- 50% lighter weight
- 140 channels

Higher SNR

- Enhanced signal-to-noise ratio
- Reduced scan time
- Improved image quality



Rotating frame of reference



Spin packet

