





Use of Magnetization Phase in Applications - 2

O. Wieben, Ph.D.
owieben@wisc.edu

Depts. of Medical Physics & Radiology
University of Wisconsin - Madison










Chemical Shift Imaging

Separate metabolites with different chemical shifts, e.g.

- Fat / Water Separation for Fat Suppression or Fat Quantification
- Metabolite maps for hyperpolarized gas imaging (^{13}C) for Lactate, Alanine, and Pyruvate









Chemical Shift Imaging


MR Spectroscopy

- High spectral information, low spatial information, long scan times
- No a priori information exploited



MR Spectroscopy



PRESS acquisition
2D multivoxel ^1H -MRSI
 32×32 voxels zero-filled
scan time: 6min53s

a: FLAIR
b: CE T1weighted

Maps:
c: Choline map
d: N-acetylaspartate
e: Creatine (Cr)
f: Lactate / Lipids

*A Di Costanzo et al.
Eur Radiol 2007*



Fig. 1 Example of 2D-multivoxel proton spectroscopic imaging (^1H -MRSI) performed using a 3.0-T MR scanner (Signa Horizon LX) with a standard head coil, and the protocol (PROBES) and the software (Functional Performance) provided by the manufacturer (General Electric Medical Systems, Milwaukee, WI). In brief, a point-resolved spectroscopy (PRESS) pulse sequence, with repetition time (TR) 1,500 ms, echo time (TE) 144 ms, field of view (FOV) 24 cm, phase 16, section thickness 10 mm and acquisition time 6 min 53 s, provided a spectroscopic data set, which was spatially zero-filled to 32×32 images with a final voxel resolution of $7.5 \times 7.5 \times 10$ mm (0.56 cm 3). Fluid-attenuated inversion recovery (FLAIR) (a) and contrast-enhanced T $_2$ -weighted (b) images from a 45-year-old man with a left fronto-temporo-parietal glioblastoma. Choline (Cho) (c), N-acetylaspartate (NAA) (d), creatine (Cr) (e) and lactate/lipids (LL) (f) maps. Localizing FLAIR image (g) and spectra (1-6) from the corresponding boxes in g. Color scale: red strongest signal intensity, blue weakest signal intensity. Higher spatial resolution improves the ability to recognize the spatial heterogeneity of high-grade gliomas. The spectrum 1, showing a low level of Cho, underdetectability of Cr and NAA, and high LL peak, suggests the presence of necrosis. The spectrum 2, in the margins of glioma, presents the typical tumor pattern, characterized by an abnormal Cho/NAA ratio (>1). The spectra 3 and 4, in an apparently edematous region surrounding the enhanced margins of glioma, show low metabolite levels, as compared to the normal spectra 5 and 6, but different pattern: the spectrum 3 has a tumor pattern and suggests the presence of infiltrating tumor cells; the spectrum 4 has a normal Cho/NAA ratio and suggests noninfiltrated vasogenic edema.


*A Di Costanzo et al.
Eur Radiol 2007*

MR-Signal Phase


Chemical Shift

- Fat-water difference = 220 Hz @ 1.5T / 440 Hz @ 3T
- Displacement fat vs. water \rightarrow **Misregistration**


1 pixel = $2 \times \text{BW} / \text{matrix Hz}$



BW = 6.94kHz
 \rightarrow 4 pixel shift



BW = 31.25kHz
 \rightarrow 0.9 pixel shift



BW = 9.23kHz
 \rightarrow 3 pixel shift

Degree of misregistration $\left\{ \begin{array}{l} \text{- Bandwidth} \\ \text{- Field strength} \end{array} \right.$

Solving for 'arbitrary' TE's

$$\hat{S}(t_n) = (\hat{\rho}_w e^{-j2\pi\Delta f_w t_n} + \hat{\rho}_f e^{-j2\pi\Delta f_f t_n}) e^{-j2\pi\Delta f t_n}$$

5 unknowns: 2×complex
1×scalar

Simple solution: Fieldmap=0 -> Matrix Inversion

Matrix Notation: $\hat{S} = A \hat{\rho}$

Better solution: iterative approach

- estimate field map ΔB_0
- no phase unwrapping required
- 'free choice' of TR und ΔTE

Fat Water Separation with IDEAL

SB Reeder, MRM, 2004

Source Water Fat

Multiecho GRE

O Wieben et al., ISMRM 2005, p. 2386

EPSI – Timing Properties

$$f_{\text{Nyquist}} = 1 / (2 \cdot \Delta TE)$$

$$\Delta f = 2 \cdot f_{\text{Nyquist}} / M = 1 / (\Delta TE \cdot M) = 1 / TE_{\text{width}}$$

$$TE_{\text{eff}} = \Delta TE \cdot (M-1)$$

EPSI – Proof of principle

Proof of principle

FLASH images (no need for ultrashort TR)

Single echo measurements (no 'EPI problems')

3 components

- Water – Acetone – Fat

Imaging Sequence

- TR = 40 ms
- 20 Echos: 5 ms ≤ TE ≤ 4 ms
- ΔT = 1 ms -> f_{Nyquist} = 500 Hz
- T_{width} = 20 ms -> Δf = 50 Hz

Daten vom 28 Jan 04 - SPGR







SPGR – Mag and Phase







Acetone Spaced

Wasser Wasser

T_{width} = 20ms, ΔTE = 1 ms
-> Δf = 50 Hz, f_{Nyquist} = 500

Daten vom 28 Jan 04 - SPGR





1-Directional Velocity Encoding

ECG gated, velocities encoded in **phase difference image $\Delta\phi$**

Magnitude

Phase Diff $\Delta\phi$

Flow, z

1-Directional Velocity Encoding

ECG gated, velocities encoded in **phase difference image $\Delta\phi$**

Magnitude

Phase Diff $\Delta\phi$

Flow, z

Desired parameters

- Single slice
- Readout: 256
- # of phase encodes: 160
- TR = 8 ms
- Views per segment = 6
- Temporal resolution = 48 ms
- Heart rate = 60 bpm
- # of cardiac phases: 20

Scan duration

$160/6 * 2 = 54$ heart beats
= 54 s
-> shorten scan time to single breath hold!
view sharing
partial Fourier
parallel imaging
...

Common Applications

Heart & Great Vessels (breath hold)

- Cardiac output
- Regurgitation volumes (valves)
- Valve insufficiency
- Stenoses
- Coarctation
- ...

Neurovascular

- Carotid arteries and others

C. Francois, University of Wisconsin

Common Applications

Heart & Great Vessels (breath hold)

- Cardiac output
- Regurgitation volumes (valves)
- Valve insufficiency
- Stenoses
- Coarctation
- ...

Neurovascular

- Carotid arteries and others

C. Francois, University of Wisconsin

3-Directional Velocity Encoding

Acquire four data sets:

- (1) with flow comp
- (2) with flow in slice
- (3) with flow in phase
- (4) with flow in read

α TE

G_{slice}

G_{phase}

G_{read}

Courtesy of Klaus Scheffler, University of Basel

Volumetric 3-directional PC MRI

Acquisition

- Volumetric coverage
- 3-directional flow encoding: 4 acquisitions
- ECG gating
- Breathing motion in abdomen and chest

Reduce acquisition times

- View sharing & advanced resp gating
- Radial undersampling (PC VIPR)
- kt BLAST
- ...

S/I A/P R/L

'4D MR Flow'

Time = 20.40 ms

Volumetric Coverage
Three-directional Velocity Encoding
Velocity Vector Field
Dynamic

L. Wigstrom et al. MRM 36(5) 1996

- Magnitude and velocity field
 -> inherently coregistered
- 10-25 min scan time
- 15-20 cardiac phases
- Spatial resolution: (1-3 mm)³
- Major advances over the last decade
 - reduced scan time

backup
 Courtesy of A. Frydrychowicz and Markl,
 University of Freiburg

4D MR Flow examples

Cranial venous flow
 - 'Streamlines'

Hepatic flow – 'Particle Traces'

K.M. Johnson et al., ISMRM 2008, p.733

Comprehensive PC MR

Vascular Anatomy

Velocity vector field
 Cardiac gating
 Volumetric Imaging

Comprehensive Information
 Vascular anatomy
 3D Velocity fields
 Hemodynamic parameters
 +noninvasive

3D Velocity Fields

Visualization

Flow measurements

Pressure gradients

Wall shear stress

Velocity Encoding

Upper limit
 - **Velocity Sensitivity (venc)**

- Max detectable flow velocity
- Velocity > venc → **aliasing**

→ Adjust venc to expected maximum velocity

venc > v_{Max}

Velocity Encoding

Upper limit
 - **Velocity Sensitivity (venc)**

- Max detectable flow velocity
- Velocity > venc → **aliasing**

→ Adjust venc to expected maximum velocity

venc > v_{Max}

small venc

Venc Settings

Magnitude

50 cm/s

30 cm/s


5 cm/s

optimal Venc for vascular flow
 too high for CSF flow

Venc too low for vascular flow
 too high for CSF flow

Venc too low for vascular flow
 optimal for CSF flow

2D phase contrast velocity encoding



Case 1. Normal Subject


Case 2. Neurologically Normal
Subject But Abnormal
Ventricle Size and CSF Flow

Case 3. Patient With Adult
Communicating Hydrocephalus



Zhu et al. J Magn Reson Imaging 24:756-770 (2006)

2D phase contrast velocity encoding



Case 1. Normal Subject

Case 2. Neurologically Normal
Subject But Abnormal
Ventricle Size and CSF Flow

Case 3. Patient With Adult
Communicating Hydrocephalus



Zhu et al. J Magn Reson Imaging 24:756-770 (2006)