


Magnetic Resonance Imaging: The RF subsystem

Madhav Venkateswaran
Krishna Kurpad
Bahareh Behzadneshad
Oliver Wieben
09/27/2016



MRI Scanner Components



MR Signal Detection


\vec{M} : Magnetization

$$\vec{M} = M_0 \vec{e}_z$$

f_0 : Larmor Frequency

$$f_0 = \frac{\gamma}{2\pi} B_0$$

$$= 63.86 \text{ MHz}$$



B_1 = Transverse Field

$$B_1(t) = \hat{B}_1 \sin(2\pi f_0 t)$$

$$\hat{B}_1 = 20 \text{ } \mu\text{T} = 0.2 \text{ G}$$

Contents

- The RF coil
 - Transmit Operation
 - Receive Operation
 - Geometry and B_1 Field Distribution
- Matching and tuning
- RF subsystem
 - Transmit
 - Receive

The radiofrequency coil

• Functions

- Set up uniform transverse RF magnetic field
- Sense NMR signal

• Types


- Surface coils
 - Single surface coils
 - Multi-channel arrays
 - Linear arrays
 - Volume arrays
- Volume coils
 - Linear
 - Quadrature



The radiofrequency coil

- Constructed using conductive wire, rods or tape
- Coil impedance is minimum at resonance




$$Z = R + j\omega L + \frac{1}{j\omega C}$$



The radiofrequency coil

- Principle of operation:
 - Transmit phase:
 - Ampere's circuital law


The integral of the magnetic field around a closed loop that encloses a current carrying conductor is proportional to the magnitude of the electric current in that conductor.

$$\oint_c B \cdot dl \propto I_{cond}$$





The radiofrequency coil

- Principle of operation:
 - Receive phase:
 - Faraday's law of electromagnetic induction


The electromotive force (emf) induced around a closed path is proportional to the rate of change of magnetic flux through the surface bounded by that path.

$$\varepsilon \propto -\frac{d}{dt} \iint_s B \cdot dS$$



Birdcage Coil – RF Excitation

$$J_n = J_s(t) \sin(2\pi(n-1)/N)$$

$$J_s(t) = \hat{J}_s \cos(2\pi f_0 t)$$


Linear Polarized Field




$$J_n = J_s(t) \sin(2\pi(n-1)/N)$$

$$J_s(t) = \hat{J}_s \cos(2\pi f_0 t)$$



$$\vec{H}(t) = \hat{H} \cos(2\pi f_0 t) \vec{e}_x$$

$$H(t) = \frac{\hat{H}}{2} (e^{j2\pi f_0 t} + e^{-j2\pi f_0 t})$$



The transmit coil...volume coil


- Axial current on infinitely long cylindrical surface whose amplitude varies as the sine of the azimuthal angle produces uniform transverse field
- Birdcage coil: approximation of sinusoidal current distribution


Single feed point:
Transverse field is linearly polarized

For circularly polarized field, second feed point placed 90° to first and excited at quadrature phase difference

The birdcage coil (Eddy Boskamp, GE Medical Systems)



Circular Polarized Field



$$J_n = J_{A,n} + J_{B,n}$$


$$J_{A,s}(t) = \hat{J}_s \cos(2\pi f_0 t)$$

$$J_{B,s}(t) = \hat{J}_s \sin(2\pi f_0 t)$$

$$\vec{H}_x(t) = \hat{H} \cos(2\pi f_0 t) \vec{e}_x$$


$$\vec{H}_y(t) = \hat{H} \sin(2\pi f_0 t) \vec{e}_y$$

$$H(t) = \hat{H} e^{j2\pi f_0 t}$$



The radiofrequency coil


- RLC circuit
- Constructed using conductive wire, rods or tape
- R and L are inherent in conductive element
- C is used to resonate the coil
- Coil impedance is minimum at resonance



Krishna N. Kuppala, 10/10/2007


Matching and Tuning

- Matching
 - PA, LNA are 50 Ω systems
 - Match coil Z to 50 Ω for max power transfer
- Tuning
 - Tune coil to desired center frequency





Most commonly used for MR coils:
2 element 'L' matching networks

C_T tunes the coil; C_M or L_M match to 50 Ω



Krishna N. Kuppala, 10/10/2007



RF coil...Match & Tune procedure

- Measure coil impedance
 - Careful about choosing the measurement plane!
- Determine series coil capacitance
 - If coil inductance is too large, reduce it to a manageable value
- Determine parallel matching capacitance
 - Impedance measurement close to the 50 Ω impedance circle
- Determine series tuning capacitance
 - Final tuning/matching by adjusting the series/parallel capacitor

Krishna N. Kuppala, 10/10/2007

The transmit coil

- High homogeneity: Uniform flip angle

$$\alpha(\vec{r}) = \gamma \int_0^{\tau} B_1(\vec{r}, \tau) d\tau$$

$$S(\vec{r}) \mu \left(\sin \{ \alpha(\vec{r}) \} \right)^3$$



$$S(\vec{r}) \mu \sin \{ \alpha(\vec{r}) \}$$

Spin Echo Sequence

Gradient Echo Sequence

- High sensitivity and Q: minimum power requirement

Krishna N. Kuppala, 10/10/2007



The receive coil



- SNR optimization*

$$SNR \sim \sqrt{\frac{Q_L}{V_{eff}}} \quad \frac{1}{Q_L} = \frac{1}{Q_0} + \frac{1}{Q_{pat}} \quad V_{eff} = \int \frac{B_1^2}{B_{1(ROI)}^2} dV$$

- Patient losses should dominate
- $Q_0 > 5Q_L$
- Coil should fit anatomy of interest
- Loss mechanisms
 - Coil losses
 - Dielectric losses
 - Radiation losses
 - RF eddy current losses

High SNR
Not very homogeneous sensitivity

* From Lecture notes on RF subsystem by Dr. Ed Boskamp
Kirsteen N. Murphey, 10/10/2007



Parallel MRI: Challenges


- Generally, Signal-to-Noise Ratio is degraded
- Hardware issues
 - Multi-channel receiver hardware
 - Multi-element coil systems
 - Dedicated image reconstruction system




90 Channel Array from MGH

RF subsystem...The receive chain

Block Diagram



Block diagram of the MR system




Steven M. Wright, Texas A&M University

- Bahar Behzadnezhad
- PhD candidate in Electrical Engineering
- Mentors:
 - Nader Behdad, Dept. of ECE
 - Alan McMillan, Dept. of Radiology



34

MRI Scanner Components



35


Proton Imaging Frequencies in MRI

Magnetic Field Strength	Frequency
0.8 T	34.06 MHz
1.5 T	63.86 MHz
3 T	127.72 MHz
7 T	298.01 MHz

36

Modes of Operation


TX, RX, TX/RX



37


Mode of Operation

Linear, Circular



38

Wave Polarization





A circularly polarized wave as a sum of two linearly polarized components 90° out of phase

39

Types of RF Coils



Volume coils (Birdcage coils, TEM coils)

	😊	😞
Birdcage Coil 	Excellent transmit magnetic field homogeneity below 3T for a large FOV	Poor SNR significant coupling between the patient and the coil at high frequencies
TEM coil 	Homogeneous magnetic field at high frequencies (7T and above)	Poor SNR

40

Types of RF Coils


Surface Coils

	😊	😞
Single element 	High SNR for receive	Poor field homogeneity for transmit for larger FOV
Array Coil 	High SNR Large FOV	Design Complications Coupling between elements

41


Single Loop RF Coil

For a target depths of interest, determined with full-wave simulations




42

B1 Field distribution




Single Element Circuit Schematic


A typical circuit schematic for a receive-only coil element



Measurements




Coil Array




Decoupling

Overlapping



- Capacitive decoupling
- Preamplifier decoupling

Specific Absorption Rate (SAR)





New Developments and Outlook

- METAMATERIALS**
 - Find artificially designed materials to tackle the high SAR levels, image inhomogeneity

- WIRELESS TECHNOLOGY**

As the amount of receivers keeps increasing, the issue of cabling is a challenge

 - High interaction between cables
 - Space
 - Weights

50

New Developments and Outlook

- Flexible and Adaptive RF Coils**
 - Close fit to the subject
 - Higher SNR
- Traveling Wave Imaging**
 - Higher field homogeneity at higher frequencies (7T and above)
- Multituned RF Coils**
 - Tuned to several frequencies
 - Detect other nuclei other than hydrogen

51